



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 8

(Os01g0656200, LOC_Os01g46760)

Product Code	CSB-MP735941OFG
Abbreviation	Os01g0656200, LOC_Os01g46760
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SN75
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSSDTSRRDH AAMAVREVL A GDRKVGTVSR SARRRRLELR RLGRTASAVA EDDAAKRVRP ASDSSSDSSE SAKVAPEPTA EVARWPACVS HGAVSVIGRR REMEDAIFVA APFLAASKEA AVEGSGVAEE EGKEEDEGFF AVYDGHGGSR VAEACRERMH VVLAEEVRVR RLLQGGGGGA DVEDEDRARW KEAMAACFTR VDGEVGGAE ADTGEQTVGS TAVVAVVGPR RIVVANCGDS RAVLSRGGVA VPLSSDHKPD RPDEMERVEA AGGRVINWNG YRILGVLATS RSIGDYLLKP YVIAEPEVTV MDRTDKDEFL ILASDGLWDV VSNDVACKIA RNCLSGRAAS KYPESVSGST AADAAALLVE LAISRGSKDN ISVVVVELRR LRSRTTASKE NGR
Source	Mammalian cell
Target Names	Os01g0656200
Protein Names	Recommended name: Probable protein phosphatase 2C 8 Short name= OsPP2C08 EC= 3.1.3.16
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.