



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 40 (Os06g0170500, LOC_Os06g07350)

Product Code	CSB-BP735942OFG
Abbreviation	Os06g0170500, LOC_Os06g07350
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SNN4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAHRLLRDAQ ADGWERSDFP IICESCLGDN PYVRMLRAEY DKECKICARP FTVFRWRPGR DARYKKTEIC QTCCCKLKNVC QVCLLDLEYG LPVQVRDTAL AINSNDAIPR SDVNREYFAE EHDKARAGI DYDSSHGKAR PNDTILKLQR TAPYYKRNRA HVCSFYVRGE CTRGAECPYR HEMPETGELS QQNIKDRYYG VNDPVALKLL GKAGEMPSLT PPDESIRTL YIGGLNNRIT EQDLRDQFYA HGEIESIRMV LQRACAFVTY TTREGAEKAA EELANKLVIK GIRLKLMWVK PQAPKPEDDE AGRQGHVAHG GMLPRAVISQ QQSGDQPQPP GMEGQQQAPS GSYYFNIPAP PGAERTLYPS MDPQRMGALV KSQEGDGKPG PQQAAQAQAS SSSGQSYMP PQYYHGQYPP YPPYGGYMP PPRMPYPPPP QYPPYQPMLA TPAQSQASSS QQPAPATLHQ AQVPPPQQT QN
Source	Baculovirus
Target Names	Os06g0170500
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 40 Short name= OsC3H40
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.