



Recombinant Mouse Protein FAM46C (Fam46c)

Product Code	CSB-MP735949MO
Abbreviation	Fam46c
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5SSF7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAEEGSSTKD SESFSVLNWD QVSRLHEVLT EVVPIHGRGN FPTLEITLKD IVQTVRGRLE EAGINVQDVR LNGSAAGHVL VKDNGLGCKD LDLIFHVALP TEAEFQLVRD VVLCSSLNFL PEGVNKLKIS PVTLKEAYVQ KLVKVCTDTD RWSLISLSNK NGRNVELKFV DSIRRFQFEFS VDSFQIILDS LLFFYDCSGN PISEHFHPTV IGESMYGDFE EAFDHLQNR IATKNPEEIR GGGLLKYSNL LVRDFRPADQ EEIKTLERYM CSRFFIDFPD ILEQQRKLET YLQNHFSDEE RSKYDYLMIL RRVVNESTVC LMGHERRQTL NLISLLALRV LAEQNIIPSA TNVTCYYQPA PYVSDGNFNN YYIAHPPITY SQPYPTWLP N
Source	Mammalian cell
Target Names	Fam46c
Protein Names	Recommended name: Protein FAM46C
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>