



Recombinant Mouse Spermatogenesis-associated protein 22 (Spata22)

Product Code	CSB-EP735953MO-B
Abbreviation	Spata22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SV06
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MKRNLNESSA RSTAGCLPVP LFNQKKRNRQ PLTSNPLQND PGVSTVSDSY GSPSFPTDWA WEAVNPEVAP LKKTVENTGQI PASASYPWRS QDSVSKSIQS NAERSQSAWG YRGNNRNTSL RTWDFRPQHK TVSPAANSEF SSCPVNLGAQ QKQFQTPEF PNLPGHKEAE VPRQTCLSKL PGSTMKGPD ASALQAFKPS FQQNPFKKTV LGDIPRENSL KEGTLHQLKE KDNSLRISA VIESMKYWRA HVQKTVLLFE ILAVLDSAVT SGPHYSKTFL MRDGKNILPC VFYEIDRELP RLIRGRVHRC VGHYDPDKNI FKCVSVRPAS ASEQKTFQAF VTIADAEMKY HTKVTNEM
Source	E.coli
Target Names	Spata22
Protein Names	Recommended name: Spermatogenesis-associated protein 22
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.