



Recombinant Human Ligand-dependent nuclear receptor-interacting factor 1 (LRIF1)

Product Code	CSB-BP735968HU
Abbreviation	LRIF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T3J3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSNNLRRVFL KPAEENSGNA SRCVSGCMYQ VVQTIGSDGK NLLQLLPIPK SSGNLIPLVQ SSVMSDALKG NTGKPVQVTF QTQISSSSTS ASVQLPIFQP ASSSNYFLTR TVDTSEKGRV TSVG TGNFSS SVSKVQSHGV KIDGLTMQTF AVPPSTQKDS SFIVVNTQSL PVTVKSPVLP SGHHLQIPAH AEVKSVPASS LPPSVQQKIL ATATTSTSGM VEASQMPTVI YVSPVNTVKN VVTKNFQNIY PKPVTEIAKP VILNTTQIPK NVATETQLKG GQHSQAAPVK WIFQDNLQPF TPSLVPVKSS NNVASKILKT FVDRKNLGDN TINMPPLSTI DPSGTRSKNM PIKDNALVMF NGKVYLLAKK GTDVLPSQID QQNSVSPDTP VRKDTLQTVS SSPVTEISRE VVNIVLAKSK SSQMETKSLS NTQLASMANL RAEKNKVEKP SPSTTNPHMN QSSNYLKQSK TLFTNPIFPV GFSTGHNAPR KVTAVIYARK GSVLQSIEKI SSSVDATTVT SQQCVFRDQE PKIHNEMAST SDKGAQGRND KKDSQGRSNK ALHLKSDAEF KKIFGLTKDL RVCLTRIPDH LTSGEGFDSF SSLVKSGTYK ETEFMVKEGE RKQQNFDDKKR KAKTNKKMDH IKKRKTENAY NAIINGEANV TGSQLLSSIL PTSDVSQHNI LTSHSKTRQE KRTEMEYYTH EKQEKGTLSN NAAEQSHFF NKNYTEDIFP VTPPELEETI RDEKIRRLKQ VLREKEAALE EMRKKMHQK
Source	Baculovirus
Target Names	LRIF1
Protein Names	Recommended name: Ligand-dependent nuclear receptor-interacting factor 1 Alternative name(s): Receptor-interacting factor 1
Expression Region	1-769
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.