



# Recombinant Human Ligand-dependent nuclear receptor-interacting factor 1 (LRIF1)

<b>Product Code</b>	CSB-EP735968HU
<b>Abbreviation</b>	LRIF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5T3J3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSNNLRRVFL KPAEENSGNA SRCVSGCMYQ VVQTIGSDGK NLLQLLPIPK SSGNLIPLVQ SSVMSDALKG NTGKPVQVTF QTQISSSSTS ASVQLPIFQP ASSSNYFLTR TVDTSEKGRV TSVGTFGNFSS SVSKVQSHGV KIDGLTMQTF AVPPSTQKDS SFIVVNTQSL PVTVKSPVLP SGHHLQIPAH AEVKSVPASS LPPSVQQKIL ATATTSTSGM VEASQMPTVI YVSPVNTVKN VVTKNFQNIY PKPVTEIAKP VILNTTQIPK NVATETQLKG GQHSQAAPVK WIFQDNLQPF TPSLVPVKSS NNVASKILKT FVDRKNLGDN TINMPPLSTI DPSGTRSKNM PIKDNALVMF NGKVYLLAKK GTDVLPSQID QQNSVSPDTP VRKDTLQTVS SSPVTEISRE VVNIVLAKSK SSQMETKSLS NTQLASMANL RAEKNKVEKP SPSTTNPHMN QSSNYLKQSK TLFTNPIFPV GFSTGHNAPR KVTAVIYARK GSVLQSIEKI SSSVDATTVT SQQCVFRDQE PKIHNEMAST SDKGAQGRND KKDSQGRSNK ALHLKSDAEF KKIFGLTKDL RVCLTRIPDH LTSGEGFDSF SSLVKSGTYK ETEFMVKEGE RKQQNFDKKR KAKTNKKMDH IKKRKTENAY NAIINGEANV TGSQLLSSIL PTSDV SQHNI LTSHSKTRQE KRTEMEYYTH EKQEKGLNS NAAEQSHFF NKNYTEDIFP VTPPELEETI RDEKIRRLKQ VLREKEAALE EMRKKMHQK
<b>Source</b>	E.coli
<b>Target Names</b>	LRIF1
<b>Protein Names</b>	Recommended name: Ligand-dependent nuclear receptor-interacting factor 1 Alternative name(s): Receptor-interacting factor 1
<b>Expression Region</b>	1-769
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.