



Recombinant Human Protein FAM102B (FAM102B)

Product Code	CSB-BP735976HU
Abbreviation	FAM102B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T8I3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMKKKKFKFK VDFELEEELSS VPFVNGVLFC KMRLLDGGSF TAESSREVVQ ANCVRWKRFK SFMCKMSASA ATGILDPCIY RVSVRKELKG GKAYAKLGFA DLNLAEFAGS GNTTRRCLLE GYDTKNTRQD NSILKVLISM QLMSGDPCFK TPPSTSMSIP IAGESESLQE DRKGGETLKV HLGADLSAK SASVPDELGA CGHSRTSSYA SQQSKVSGYS TCHSRSSSFS ELCHRRNTSV GSTSTGVESI LEPCDEIEQK IAEPNLDTAD KEDTASEKLS RCPVKQDSVE SQLKRVD DTR VDADDIVEKI LQSQDFSLDS SAE E E E G L R L F V G P G G S T T F G S H H L P N R V G S GAYEQVVIKR
Source	Baculovirus
Target Names	FAM102B
Protein Names	Recommended name: Protein FAM102B
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.