



Recombinant Dictyostelium discoideum Protein VASP homolog (vasp)

Product Code	CSB-EP735989DKK-B
Abbreviation	vasp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TJ65
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MSETAIFNAT GQVFTYSPQT RNWVPSSNVP ATLQMYFNLSG ANTYRVIGRA GDDPNNFLIN FAVKSEVVYS RASEIFHQFT DQRTHFGINF TSKQDADTFG GGFENVLRSL KGGPQQPPPQ VPKPQQQPP QPQQPPQRPP STVIKPVAP QAPVAPPQAP AAPAQAPAPP AAPPAPPKPP GPPPPPPAPK PPAAGGGTGR NALLGSIENF SKGGLKKTVT VDKSAGVPTK TTPSANSAS NAGSEPSSGG STPAPAPKSS GGGGGDLMA EVMAKRAKMK AAASQPKEES SAPTPAPTPA PTPAPTPSST SFKPPQSFSK PATKTAAANK PPSPSLSAPL PSTVANEDLQ SLKEEILTEV RKEIQKAKDE ILEAIRASQH
Source	E.coli
Target Names	vasp
Protein Names	Recommended name: Protein VASP homolog Alternative name(s): DdVASP
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.