



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R735 (MIMI_R735)

Product Code	CSB-MP736041ADAZ
Abbreviation	MIMI_R735
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UNY9
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MSSEPVTYYM IVNKS WLKII NYESP KSKRS LSGY LIKLQS NHCIE KRGYF FTDKV HLLDF IDEGDTIVSI KLPQN FDEL VVRDG KGWRS NQII ENT SK LHDINTAKKF HIVGINDQQV ANLAIKYIDV SILYYLISTG VELKSNIGYI HWWNKICMDF NNPERYFVRK LLDLTPFYVV SHGPSRNKNS QTPIRAPFET QVHNPISGQE CYVQNFYFTD ADRAMKFLSN GLSVFVYKIK PYLSNNDGIN DIIREVRPNI FWATCVVNG IYDLAKEDSY DSLLDLEFNI EDTIFALKK NWAKIVYLLV SRHNFDTIFK TKLLFTAAAL GYENMIYLLM KHQPIYSASQ IIMGSELSSS FIMHHELDIY KVADLPVPER FIDTFLVLAT VGG SIRT LTA FIKKYTEHFH KHKQIMIMYA KYYSKPDLLA LFEYYSNNIV DNHITQ
Source	Mammalian cell
Target Names	MIMI_R735
Protein Names	Recommended name: Uncharacterized protein R735
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.