



# Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R735 (MIMI\_R735)

<b>Product Code</b>	CSB-YP736041ADAZ
<b>Abbreviation</b>	MIMI_R735
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UNY9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acanthamoeba polyphaga mimivirus (APMV)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSEPVTYYM IVNKS WLKII NYESP KSKRS LSGYLIK LQS NHCIEKRGYF FTDKVHLLDF IDEGDTIVSI KLPQNFDELV VVRDGKGWRS NQIIIENTSK LHDINTAKKF HIVGINDQQV ANLAIKYIDV SILYYLISTG VELKSNIGYI HWWNKICMDF NNPERYFVRK LLDLTPFYVV SHGPSRNKNS QTPIRAPFET QVHNPISGQE CYVQNFYFTD ADRAMKFLSN GLSVFVYKIK PYLSNNDGIN DIIREVRPNI FWATCVVNG IYDLAKEDSY DSLLDLEFNI EDLTIFALKK NWAIVYLLV SRHNFDTIFK TKLLFTAAAL GYENMIYLLM KHQPIYSASQ IIMGSELSSS FIMHHELDIY KVADLPVPER FIDTFLVLAT VGG SIRT LTA FIKKYTEHFH KHKQIMIMYA KYYSKPDLLA LFEYYSNNIV DNHITQ
<b>Source</b>	Yeast
<b>Target Names</b>	MIMI_R735
<b>Protein Names</b>	Recommended name: Uncharacterized protein R735
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.