



Recombinant Acanthamoeba polyphaga mimivirus

Putative ankyrin repeat protein R598 (MIMI_R598)

Product Code	CSB-MP736051ADAZ
Abbreviation	MIMI_R598
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UP65
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDDQPINIVT DIKNAICKDD LDSFIILMEN NPSINLFCEP DKQSIDFGGI IQPIEFCKNN DLIKMVNEVD DQLIDHYIST VLREEVIKNC SVKILKYLF MGLPVDFGRN FSIRLCSTKY LGSNHYIKEN PGQDTLSLLR LLIEYGVDVN AHNYLPLYSA VSSKNFDKVK LLVENGANVL RVANGNENSF YFKIDSIKYL LDNGVEIDMN LSRALFLSIK DNSIECIQFY LELGADINKI SPMDIIGLIK LRNIVGIRLL IENGFDFTSL NKLAGNNNER IIDVLINEAN VDISALLKIL LEFW
Source	Mammalian cell
Target Names	MIMI_R598
Protein Names	Recommended name: Putative ankyrin repeat protein R598
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.