



# Recombinant Acanthamoeba polyphaga mimivirus Putative TPR repeat-containing protein R856 (MIMI\_R856)

<b>Product Code</b>	CSB-BP736120ADAZ
<b>Abbreviation</b>	MIMI_R856
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UQQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acanthamoeba polyphaga mimivirus (APMV)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEIINACLKN VGVEENMIEK VIETFQKEND ECYNLVHIFN KAAIVFHRNG QHKSLEMYE KAFGNIFNGE FALSDFYSV NGMASMYQAL GDYDIAIKKY NSVIKIKDM CLDNNSDLVY ALMGIASISQ IKGNYDEALS KYNEALEINE KLYGRNHiet AFVLNRLGML YHELDDNDKS IDHFNESLKI YREKYPNKLf NIAFTISRLA QSLKMGNDs EALEKYQESi DIFNKIFTIS HQAVAFSLYG IGTVYEFrSE YSKALEKYQE SLQTYKNVYe RSEKYQHYSi ASCLYKIGLV YKLSGNDNES TTYLNQANQM FESTSTNIND KNYQACKKFL QD
<b>Source</b>	Baculovirus
<b>Target Names</b>	MIMI_R856
<b>Protein Names</b>	Recommended name: Putative TPR repeat-containing protein R856
<b>Expression Region</b>	1-342
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.