



# Recombinant Haloarcula marismortui Signal recognition 54 kDa protein (srp54)

<b>Product Code</b>	CSB-MP736158HTJ
<b>Abbreviation</b>	srp54
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UY20
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Haloarcula marismortui (strain ATCC 43049 / DSM 3752 / JCM 8966 / VKM B-1809) (Halobacterium marismortui)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVLDDLGSSL RGTLDDLRGK SRLSEEDIED IVKEIQRSLQ QADVVDVGLVQ DLSNSIETRA LDEEPPAGTT PRDWVLRIVY EELVDLVGES TELPLEEQTI MLAGLYGSGK TTTAAKMAWW FSTKGLRPAI IQTDTDRPGA YDQSKEMAER AEVDFYGDPE EDDPVKIARD GLEATENADV RIVDTAGRDG LNEELIEQIE RIEQEVQPDR DLLVLDAAMG QSAKSQAADF EAAIGIDGVV ITKLDGTAKG GGALAAVNET DSTIAFLGSG ETVKDIERFE PSGFISRLLG MGDCLKQLTER VERAMEETQE GDEEDWDPEM MLEGQFTLKD MRKQMQTMNN MGPLDQVMDM IPGLGGGLMD QLPDDAMDVT QERMQDFDVI MDSMTEEELE NPRVVGQSRT KRICRGSQKP EERVRELLQQ HKQMEQMLKQ FQGMGDGME RMMKMQQGG GGGGGMGGMG GGGMGPFGD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	srp54
<b>Protein Names</b>	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
<b>Expression Region</b>	1-469
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.