



Recombinant Haloarcula marismortui CCA-adding enzyme (cca)

Product Code	CSB-EP736166HTJ-B
Abbreviation	cca
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5V030
Product Type	Recombinant Protein
Immunogen Species	Haloarcula marismortui (strain ATCC 43049 / DSM 3752 / JCM 8966 / VKM B-1809) (Halobacterium marismortui)
Purity	≥85% (SDS-PAGE)
Sequence	MSDEFDAVVG KVRARASPTD DERAQLQRVA DAVMADAEAA IADLPVEAEV VQVGSTARGT WTAGDRDVDV FVCFPPSIDR EALEEYGLAV GHDVLPDGRE EYAEHPYVVG EREGYAVDLV PCYAVENATE IQSAVDRTPF HTRYLQERLD DNSAAEVRVA KQFLKGIGVY GSDLRTRGFS GYLTELLVLE FGGFRAFLEA VADWHPPVRL DPDDHGSETF DDPLVVIDPT DPERNVA AVL SETNVATLQH YARDLLAEPR VSLFTEDDPC PFEAADVEAA VSQRGTTTPVA LRFAAPDVVD DQLWPQLRKS LDGLCSELDR RGFEVLRSA FVEDDSGKPE TLDTESRGRD VLLLEFAVA EQPAVERHEG PPVHVREHAS GFFQKYDDNS EVAGPFIDGD RYVVERQRAF TTATGFLSSA AVYDVGLGQR IESALENGYE VLVGTDIAAL ADGFGVDLAS YFDPKP
Source	E.coli
Target Names	cca
Protein Names	Recommended name: CCA-adding enzyme EC= 2.7.7.72 Alternative name(s): CCA tRNA nucleotidyltransferase tRNA CCA-pyrophosphorylase tRNA adenylyl-/cytidyl- transferase tRNA nucleotidyltransferase tRNA-NT
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.