



Recombinant Subtilisin-like protease 6 (SUB6)

Product Code	CSB-EP736187TBAV-B
Abbreviation	SUB6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VJ72
Product Type	Recombinant Protein
Immunogen Species	Trichophyton verrucosum (Cattle ringworm fungus)
Purity	≥85% (SDS-PAGE)
Sequence	THQ DNVPSWGLAR VSTRKPGGTT YYYDPSAGKG VTAYVIDTGI DTKHEDFGGR AKWGKNLVDQ MDEDCNGHGT HVAGTVGGTK YGLAKGVSLV AVKVLDCESG GSNSGVIKGM EWAMMDASGG GNGTAKAAGK AVMNMSLGGP RSEATNQAAK AISDAGIFLA VAAGNENMDA QHSSPASEPS VCTVAASKD DGKANFSNFG SVVDVYAPGK DIVSLKPGGG TDTLSGTSMASPHVCGLGAY LIGLGKQGGP GLCDTIKEMA NDAIQSPGED TTSKLIYNGS GK
Source	E.coli
Target Names	SUB6
Protein Names	Recommended name: Subtilisin-like protease 6 EC= 3.4.21.-
Expression Region	128-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.