



Recombinant Human CUB domain-containing protein 2 (CDCP2)

Product Code	CSB-EP736205HU-B
Abbreviation	CDCP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VXM1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QAMEGVKC GGVLSAPSGN FSSPNFPRLY PYNTECSWLI VVAEGSSVLL TFHAFDLEYH DTCSFDFLEI YNGASPDKGN LLGRFCGKVP PPPFTSSWHV MSVIFHSDKH VASHGFSAGY QKDVC GGVL T GLSGVLTSPE YPNNYPNSME CHWVIRAAGP AHVKLVFVDF QVEGNEECTY DYVAVLGGPG PTRGHHYCGS TRPPTLVSLG HELQVVKSD FNIGGRGFKA YYFSGECQEV YMAMRGNFSS PQYPSSYPNN IRCHWTIRLP PGYQVKVFFL DLDLEEPSL TKTCDFDHLA AFDGASEEAP LLGNWCGHHL PPPVTSSHNQ LLLLLHTDRS TTRRGFSVAY IGGQLGCGSG STEGEGEALQ PQSLQSPSSI PPVCPAPPMN GLLQLLHLWL HPCPLSGPLR LDGTAPACFH YCRASFPSF
Source	E.coli
Target Names	CDCP2
Protein Names	Recommended name: CUB domain-containing protein 2
Expression Region	23-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.