



Recombinant Human MORN repeat-containing protein 5 (MORN5)

Product Code	CSB-YP736208HU
Abbreviation	MORN5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VZ52
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEYTGSKYIG EYVDGRMEGK AKYILPTETI YVGEMKDGMF HEGTLYFPS GSQYDAIWEN GLAIKGTYTF SDGLHYDEKN WHYCDGYDRR FYTEILNGLK PAGMAQLTNM DPPRKIPKGY YDCGDGFYNP VTRVVKDYRN RFLRNADDDE HEWITRTRCK G
Source	Yeast
Target Names	MORN5
Protein Names	Recommended name: MORN repeat-containing protein 5
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.