



# Recombinant *Bacillus clausii* Ferrochelatase (hemH)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP736265BAAC                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Abbreviation</b>      | hemH                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q5WHT1                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | <i>Bacillus clausii</i> (strain KSM-K16)                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MEKKKIGLLV MAYGTPRKKE EIEPYYYTHIR HGRKPSEEAL ADLQERYEAI<br>GGTSPLARIT DEQANGLAAV LNEQQEDIVF VPYLGLKHID PFIEDAVQKM<br>KNDGIETAVS IVLAPHFSTF SVQSYNKRAQ EEADRIGGPT IYTIESWYKE<br>PGFLQFWKEQ IENEMKKADS IDDLCVIFSA HSLPEKIIQA NDPYPEQLRE<br>TADLLADMAN VPNYEIGWQS EGNTPEPWLG PDVQDLTREL YAKKGYKTFM<br>YCPVGFAEH LEVLYDNDYE CKLVTDELGV RYLRPMPNA DTRFLKALAD<br>VVTKRLAAVV TK                                    |
| <b>Source</b>            | <i>E.coli</i>                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Target Names</b>      | hemH                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: Ferrochelatase EC= 4.99.1.1 Alternative name(s): Heme synthase Protoheme ferro-lyase                                                                                                                                                                                                                                                                                                      |
| <b>Expression Region</b> | 1-312                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |