



Recombinant Legionella pneumophila Orotidine 5'-phosphate decarboxylase (pyrF)

Product Code	CSB-EP736344LAAB-B
Abbreviation	pyrF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5X5E0
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Paris)
Purity	≥85% (SDS-PAGE)
Sequence	MTPKLIVALD FDNQDNALQL VDKLDPNHCA LKVGSELFTEL LGPQFVKELV RREFKVFLDL KFHDIPNTVA KACHSAAELG VWMINVHAIG GLRMLQAARE SLKTYGKDRP LLIATVLTLS FEEGELASVG ISNTLPEQAT HLAMLAREAG LDGVVSSAHE VKIHKQKCGE NFITVTPGIR LPNNLKDDQS RVMTPQQAIR EGSDFLVIGR PITQASNPHE VVSALLRDL
Source	E.coli
Target Names	pyrF
Protein Names	Recommended name: Orotidine 5'-phosphate decarboxylase EC= 4.1.1.23 Alternative name(s): OMP decarboxylase Short name= OMPDCase Short name= OMPdecase
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.