



Recombinant 6-phosphofructokinase (pfkA)

Product Code	CSB-BP736393SMZ
Abbreviation	pfkA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5XBV2
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MKRIAVLTSG GDAPGMNAAI RAVVRKAISE GMEVYGINRG YAGMVDGDIF PLGSKEVGDK ISRGGTFLYS ARYPEFAQLE GQLAGIEQLK KHGIEGVVVI GGDGSYHGAM RLTEHGFPVAV GIPGTIDNDI AGTDYTIGFD TAVNTAVEAI DKLRDTSSSH GRTFVVEVMG RNAGDIALWA GIASGADQII VPEEEFDIEK VASTIQYDFE HKGKNHHIIV LAEGVMSGEA FAQKLKEAGD KSDLRVTNLG HILRGGSPATA RDRVIASWVG SHAVELLKEG KGGLAVGIHN EELVESPILG TAEEGALFSL TEEGKIIVNN PHKARLDFAA LNRSLSQ
Source	Baculovirus
Target Names	pfkA
Protein Names	Recommended name: 6-phosphofructokinase Short name= Phosphofructokinase EC= 2.7.1.11 Alternative name(s): Phosphohexokinase
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>