



Recombinant Peptide chain release factor 1 (prfA)

Product Code	CSB-MP736396SMZ
Abbreviation	prfA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XC69
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MNIYDQLQAV EDRYEELGEL LSDPDVVSdT KRFMELSREE ANTRETVTAY REYKQVIQTI SDAEEMIKDA SGDPELEEMA KEELKESKAA KEEYEEKLKI LLLPKDPNDD KNIILEIRGA AGGDEAALFA GDLLTMYQKY AETQGWRFEV MESSVNGVGG IKEVVAMVSG QSVYSKLYE SGAHRVQRVP VTESQGRVHT STATVLVMP E VEEVEYDIDP KDLRVDIYHA SGAGGQNVNK VATAVRMVHI PTGIKVEMQE ERTQQKNRDK AMKIIRARVA DHFAQIAQDE QDAERKSTVG TGDRSERIRT YNFPQNRVTD HRIGLTLQKL DTILSGKMNE VIDALVMYDQ TKKLESLNN
Source	Mammalian cell
Target Names	prfA
Protein Names	Recommended name: Peptide chain release factor 1 Short name= RF-1
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.