



Recombinant Enolase (eno)

Product Code	CSB-BP736407SMZ
Abbreviation	eno
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5XD01
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	SIITDVYAR EVLDSRGNPT LEVEVYTESG AFGRGMVPSG ASTGEHEAVE LRDGDKSRYL GLGTQKAVDN VNNIIAKAII GYDVRDQQAI DRAMIALDGT PNKGKLGANA ILGVSIAVAR AAADYLEVPL YTYLGGFNTK VLPTPMMNII NGGSHSDAPI AFQEFMIMPV GAPTFKEGLR WGAEVFHALK KILKERGLVT AVGDEGGFAP KFEGTEDGVE TILKAIEAAG YEAGENGIMI GFDCASSEFY DKERKVYDYT KFEGEGAAVR TSAEQVDYLE ELV NKYPIIT IEDGMDENDW DGWKVLTERL GKR VQLVGDD FFVTNTEYLA RGIKENAANS ILIKVNQIGT LTETFEAIEM AKEAGYTAVV SHRSGETEDS TIADIAVATN AGQIKTGSLS RTDRIAKYNQ LLRIEDQLGE VAQYKGIKSF YNLKK
Source	Baculovirus
Target Names	eno
Protein Names	Recommended name: Enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2-phosphoglycerate dehydratase
Expression Region	2-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>