



Recombinant Probable dipeptidase A (pepDA)

Product Code	CSB-MP736408SMZ
Abbreviation	pepDA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XD16
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MACTTILVGK KASYDGSTMV ARTEDSQNGD FTPKKMIVVK PEDQPRHYRS VQSSFEMDLP DNPMTYTSVP DALGKDCIWA EAGVNEANVT MSATETITTN SRVLGADPLV ASGIGEEDMV TLVLPYIRSA REGVLRGAI LEDYGTYESN GVAFSDEHDI WWLETIGGHH WIARRVPDDA YVTNPNQFGI DHFEFNNPED YLCSADLKDF IDTYHLDLTY SHEHFNPXYA FGSQRDKDRH YNTPRAWIMQ KFLNPEIVQD PRSFALAWCQ KPYRKITVED VKYVLSSHYQ DMGYDPYGSE GTPVSKKVFR PIGINRTSQT AILHIRPNKP QEIAAIQWMA YGSMFPNTMV PFFTQVKTIP DYFANTYENV STDNFYWTNR LIAALADPHY NHHETDLNRY LEETMAKGHA MLHAVEAQLL AGETVDLEEE NQKMSDYVQG ETQTLLNKIL FDASNLMTNR FSLSD
Source	Mammalian cell
Target Names	pepDA
Protein Names	Recommended name: Probable dipeptidase A EC= 3.4.-.-
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.