



Recombinant Adenylosuccinate synthetase (purA)

Product Code	CSB-YP736419SMZ
Abbreviation	purA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XE45
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MTSVVVVGTTQ WGDEGKGKIT DFLSADAEVI ARYQGGDNAG HTIVIDGKKF KLHLIPSGIF FPQKISVIGN GVVVNPKSLV KELAYLHDEG VTTDNLRISD RAHVILPYHI QLDQLQEDAK GDNKIGTTIK GIGPAYMDKA ARVGIRIADL LDKDIFAERL RINLAEKNRL FEKMYDSTPL DFDAIFEEYY AYGQEIKQYV TDTSVILNDA LDAGKRVLFE GAQGVMLDID QGTYPFVTSS NPVAGGVTTIG SGVGPSKINK VVGVCAYTS RVGDGPFPTTE LFDEVGERIR EVGHEYGTTT GRPRRVGWFD SVVMRHSRRV SGITNLSLNS IDVLSGLDAV KICVAYDLTG ERIDYYPASL EQLKRCKPIY EELPGWQEDI TGVRSLELDP ENARNYVRRRI GELVGVRIST FSVGPGREQT NILES VWASI
Source	Yeast
Target Names	purA
Protein Names	Recommended name: Adenylosuccinate synthetase Short name= AMPSase Short name= AdSS EC= 6.3.4.4 Alternative name(s): IMP--aspartate ligase
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.