



Recombinant Phosphate acyltransferase (plsX)

Product Code	CSB-BP736424SMZ
Abbreviation	plsX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XEF9
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MKRIAIDAMG GDNAPKAIVE GVNQVIEAFS DIEIQLYGDQ TKINSYLIQS DRVAIIHTDE KIMSDDEPAK AVRRKKKASM VLAAKAVKEG KADAIISAGN TGALLAVGLF VVGRIKGVDR PGLLSTIPTV TGLGFDMLDL GANAENTVKH LHQYAILGSF YAKNVRGIAN PRVGLLNNGT EETKGDPLRK ATYELLTADN TISFVGNVEA RELMSGVADV IVSDGFTGNA VLKSIEGTGI SIMGQLKQII NSGGIKTKIG ASLLKSSLYE MKKTLDYSSA GGAVLFGLKA PVVKSHGSSD VKAIFSTIKQ VRTMLDTNVV GQLVEEFAKE TQVND
Source	Baculovirus
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]--phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.