



Recombinant Rat Ribonucleoprotein PTB-binding 1 (Raver1)

Product Code	CSB-MP736445RA
Abbreviation	Raver1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XI28
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	AADVSVTHR PPLSPEAEAE AETPETVDRR TPEQELPPLD PEEIRKRLEH TERQFRNRK ILIRGLPGDV TNQEVHDLLS DYELKYCFVD KYKGTAFVTL LNGEQAEAAI NTFHQ SRLRE RELSVQLQPT DALLCVANLP PSLTQAQFEE LVRPFGSLER CFLVYSERTG HSKGYGFAEY MKKDSAARAK SDLLGKPLGP RTLYVHWTDA GQLTPALLHS RCLCVDHLPP GFNDVDALRQ ALSAVYTPTF CQLASGQDGQ LKGFVLEYE TAEMAEAAQQ RADGLALGGS HLRVSFCAPG PPGRSMLAAL IAAQATALNR GKGLLEPNI LQLLNNLGPS ASLQLLLNPL LHGGASGKQG LLGAPPAMPL LSGPALSTAL LQLALQSQNQ SQQSQQKKPG ILGDSPLGTL QAGAQPNSL LGELSAGGGL APELPPRRGK PQPLLPPLL PSGGDREPMG LGPPASQLTP PPAPMGLRGS SLRGLPKDSG PLPTPPGVSL LGEPKDYRI PLNPYLNLS LLPSSNLGK ETRGWGSSGR GRRPAEPLP SPAVPGGNA SNNGSKAFPM KPRLLSPIAS NRLPPEGLP DSYSFDYPTD VGPRRLFSHP RESNLGAHGP SRHKMSPPPS SFSEPRSGGG SGGPLSHFYS GSPTS YFTSG LQAGLKQSHL NKA VGSSPMG SSEGLLGLGP GPNGHSHLLK TPLGGQKRSF SHLLPSPEPS PEGSYVGQHS QGLGGHYADS YLKRKRIF
Source	Mammalian cell
Target Names	Raver1
Protein Names	Recommended name: Ribonucleoprotein PTB-binding 1 Alternative name(s): Protein raver-1
Expression Region	2-748
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.