



Recombinant Rat CCR4-NOT transcription complex subunit 10 (Cnot10)

Product Code	CSB-MP736452RA
Abbreviation	Cnot10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XIA4
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	AADKPADQG AEKHEGAGQS SGVTDQEKEL STSAFQAFTS GNYDACLQHL ACLDINKDD YKIILNTAVA EFFKNNQTTT DNLRQTLNQL KNQVHSAVEE MDGLDDVENS MLYYNQAVIL YHLRQHTEAI AVGEKLYQFI EPFEEKFAQA VCFLLDLYI LTHQAEKALH LLAVLEKMIS QGSGNKNGKS ETGNNSSKDG SNPKAESAL IEAAKSKIHQ YKVRGYIQMK SLKACKREIK SVMNTAGNSA PSLFLKSNFE YLRGNYRKAV KLLNSSNIAE HPGFMKTGEC LRCMFWNNLG CIHFAMSKHN LGIFYFKKAL QENDNVCAQL SAGNTDPGKK FSGRPMCTLL ANKRYELLYN CGIQLLHVGR PLAAFECLVE AVQVYHANPR LWLRLAECCI AANKGTSEQE TKGLPTKKGI VQSIVGQGYH RKIVLASQSI QNTVYNDGQS SAIPVASVEF AAICLRNALL LLPEEQQDPK QENGSKSSSQ LGGNAESSES SETCSSKSHD GDKLIPAPPS SPLRKQELN LKCSILACSA YVALALGDNL MALNHADQLL QQPKLSGSLK FLGHLYAAEA LISLDRISDA ITHLNPENVT DVSLGISSNE QDQGSCKGEN EAMESSGKRA PQCYPSSVNS ARTVMLFNLG SAYCLRSEYD KARKCLHQAA SMIHPKEVPP EAILLAVYLE LQNGNTQLAL QMIKRNQLLP AVKAHSDVRK KTVFQPVHPI QPIQMPAFTT VQRK
Source	Mammalian cell
Target Names	Cnot10
Protein Names	Recommended name: CCR4-NOT transcription complex subunit 10
Expression Region	2-744
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.