



Recombinant Rat NADH dehydrogenase [ubiquinone] iron-sulfur protein 4, mitochondrial (Ndufs4)

Product Code	CSB-BP736455RA
Abbreviation	Ndufs4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XIF3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	ADGQTRDT QLITVDEKLD VTPLTGVPEE HIKTRKVRIF VPARNNMQSG VNNTKKWKME FDTREWRWENP LMGWASTADP LSNMVLTFSA KEDAVAFSAEK HGWSYDVEGR KVPKPKSKSY GANFSWNKRT RVSTK
Source	Baculovirus
Target Names	Ndufs4
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 4, mitochondrial Alternative name(s): Complex I-18 kDa Short name= CI-18 kDa NADH-ubiquinone oxidoreductase 18 kDa subunit
Expression Region	43-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes an accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), or NADH:ubiquinone oxidoreductase, the first multi-subunit enzyme complex of the mitochondrial respiratory chain. Complex I plays a vital role in cellular ATP production, the primary source of energy for many crucial processes in living cells. It removes electrons from NADH and passes them by a series of different protein-coupled redox centers to the electron acceptor ubiquinone. In well-coupled mitochondria, the electron flux leads to ATP generation via the building of a proton gradient across the inner membrane. Complex I is composed of at least 41 subunits, of which 7 are encoded by the mitochondrial genome and the remainder by nuclear genes.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.