



Recombinant *Oryza sativa* subsp. *japonica* Zinc finger CCCH domain-containing protein 46 (Os06g0704300, LOC_Os06g49080)

Product Code	CSB-BP736540OFG
Abbreviation	Os06g0704300, LOC_Os06g49080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5Z807
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MSRRQEICRN FQRGSCKYGA QCRYLHASPH QQQQQQQAQK NPFGFGTGSR QQQQPSFGSQ FQQQQQQQK PNPFGFGVQG ANAQSRNAPG PAKPFQNKWV RDPSAPTKQT EAVQPPQAQA AHTSCEDPQS CRQQISEDFK NEAPIWKLTG YAHLRNGPCN IKGDISFEEL RAKAYEEGKQ GHSLQSIVEG ERNLQNAKLM EFTNLLNSAR PSQTPSFPTM SSFPEVKNN SFGASQTNGP PVFSSFSQIG AATNIGPGPG TTAPGMPASS PFGHPSSAPL APTFGSSQM KFGVSSVFGN QGSGQPFGSF QAPRFPSSKS PASSVQHRDI DRQSQELLNG MVTPPSVMFE ESGVNNKNEN QDDSIWLKEK WAIGEIPLDE PPQRHVSHVF
Source	Baculovirus
Target Names	LIC
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 46 Short name= OsC3H46
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.