



Recombinant Chicken Endophilin-B1 (SH3GLB1)

Product Code	CSB-EP736560CH-B
Abbreviation	SH3GLB1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5ZIR1
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MNIMDFNMKK LAADAGTFLS RAVQFTEEKL GQAEKTELDA HLENLLSKAE CTKQWTEKIM KQTEVLLQPN PNARIEEFFY EKLDRKAPSR MNNPELLGQY MIDAGNEFGP GTAYGNALIK CGETQKQIGT ADRELIQ TSA INFLTPLRNF IEGDYKTITK ERKLLQNKRL DLDAAKTRLK KAKVAEAAA SEQEVKITQS EFDRQAEITR LLEGISSTH AHHLRCLNDF VEAQMTYYAQ CYKYMLDLQK QLGSFPSTFL SNNNQSSSAP VQSVSTPSVL ASASASLPSV SNSVVTSGIS ELKSSSGSRK ARVLYDYDAA NSSELSLLAD EVITVYSIPG MDSDWLMGER GNQKKGKVPIT YLELLN</pre>
Source	E.coli
Target Names	SH3GLB1
Protein Names	Recommended name: Endophilin-B1 Alternative name(s): SH3 domain-containing GRB2-like protein B1
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>