



# Recombinant Chicken Peroxiredoxin-6 (PRDX6)

<b>Product Code</b>	CSB-EP736570CH
<b>Abbreviation</b>	PRDX6
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q5ZJF4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	PGLLLGDEA PNFEADTTQG GIRFHDFLGD SWGILFSHPR DFTPVCTTEL GRAAKLAPEF SKRNVKMIAL SIDSVPDHLA WSKDINAYNG DQPVEKLPFP IIADKDRELA VKLGMLDPDE RDKDGMPLTA RVVFIFGPKD KKLKLSILYPA TTGRNFDEIL RVVDSLQLTA YKKVATPVDW KCGDSVMVVP TLPDEEAKKL FPKGVTKDL PSGKKYLRYT PQPE
<b>Source</b>	E.coli
<b>Target Names</b>	PRDX6
<b>Protein Names</b>	Recommended name: Peroxiredoxin-6 EC= 1.11.1.15
<b>Expression Region</b>	2-224
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	<p>This protein is a member of the thiol-specific antioxidant protein family. This protein is a bifunctional enzyme with two distinct active sites. It is involved in redox regulation of the cell; it can reduce H<sub>2</sub>O<sub>2</sub> and short chain organic, fatty acid, and phospholipid hydroperoxides. It may play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury.</p>
<b>Reconstitution</b>	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>