



# Recombinant Chicken Serine/threonine-protein kinase 4 (STK4)

<b>Product Code</b>	CSB-EP736571CH
<b>Abbreviation</b>	STK4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5ZJK4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METVQLRNPR RQLKKLDEDS LTKQPEEVFD VLEKLGEGSY GSVFKAIHKE TGQVVAIKQV PVESDLQEII KEISIMQQCD SPHVVKYYGS YFKNTDLWIV MEYCGAGSVS DIIRLRNKTLL TEEIATIVQ STLKGLEYLH FMRKIHRDIK AGNILLNTEG HAKLADFGVA GQLTDTMAKR NTVIGTPFWM APGVIQEIGY NCVADIWSLG ITAIEMAEGK PPYADIHMPR AIFMIPTNPP PTFRKPELWS DAFTDFVKQC LVKSPEQRAT AIQLLQHPFV KSAKGASILR DLINEAMDIK LKRQEAQQRE LDQEDEENSE EDETDSGTMV RASGDETGTI RVVNTMSDGA NTMIEHDGTL ESQLGTMVIN TEDEEEEGTM KRRDETMQPA RPSFLEYFEQ KAKENQVNSF GKNVSGQTKN SSDWKVPQDG DYEFLKTWSV DELQRRLSAL DPMMEQEIEE IRQKYQSKRQ PILDAIEAKK RRQQNF
<b>Source</b>	E.coli
<b>Target Names</b>	STK4
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase 4 EC= 2.7.11.1 Cleaved into the following 2 chains: 1. Serine/threonine-protein kinase 4 37kDa subunit Short name= 2. MST1/N 3. Serine/threonine-protein kinase 4 18kDa subunit
<b>Expression Region</b>	1-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a cytoplasmic kinase that is structurally similar to the yeast Ste20p kinase, which acts upstream of the stress-induced mitogen-activated protein kinase cascade. The encoded protein can phosphorylate myelin basic protein and undergoes autophosphorylation. A caspase-cleaved fragment of the encoded protein has been shown to be capable of phosphorylating histone H2B. The particular phosphorylation catalyzed by this protein has been correlated with apoptosis, and it is possible that this protein induces the chromatin



condensation observed in this process.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.