



Recombinant Chicken High mobility group protein 20A (HMG20A)

Product Code	CSB-MP736583CH
Abbreviation	HMG20A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZKF4
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MENLVAGSTL PPLFADEDGS KEGSDVTVTG LAHSESSFTG GTSQPVNNPD LVEDLSQVQQ LQNESTNTAE NTEQKPEEEQ QRTKRGGWAK GRKRKKPLRD SNAPKSPLTG YVRFMNERRE QLRAKRPEVP FPEITRMLGN EWSKLPPEEK RRYLDEADRD KERYMRELEQ YQKTEAYKVF SRKAQDRQKG KLHRQD GARQ PVHDHEKEAD TKERSVFDIP IFTEEFLNHS KAREAELRQL RKS NMEFEER NAALQKHVES MRTAVEKLEV DVIQERSRNT VLQQHLET LR QALTTSFAGV PLPGSGETPT METIDSYMN R LHGIIMANPQ ENENLIATVR DVVNR LER
Source	Mammalian cell
Target Names	HMG20A
Protein Names	Recommended name: High mobility group protein 20A Alternative name(s): HMG box-containing protein 20A
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.