



Recombinant Chicken Basic leucine zipper and W2 domain-containing protein 2 (BZW2)

Product Code	CSB-MP736595CH
Abbreviation	BZW2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZL42
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MNKNQKPVLT GQRFKTRKRD EKEKFEPTVF RDTIVQGLNE AGSDLEIAIK FLDSAGSRLD YRRYADTLFD ILVAGSMLAP GGTRIDDNDK TKMTRHCVFF AEEDHDAIRN YAQVFNKLIR RYKYLEKAFE DEIKKLLLFL KAFSETEQTK LAMLSGILLA NGTLPATILT SLFTDNIVKE GIAASFAVKL FKAWMAEKDA NSVTSALRKA NLDKRLLLELF PANRQNVDFH AKYFTEAGLK ELSDFLRVQQ SLGTRKELQK ELQERLSQEC PIKEMVLYVK EEMKRNELPE PAVIGLLWTC VMNAVEWNKK EELVAEQALK HLKQYAPLLA VFSTQGQSEL ILLQKVQEYC YDNIHFMKAF QKIVVLFYKA DVLSEEAILK WYKEAHVAKG KSVFLDQMKK FVEWLQNAEE EFRI
Source	Mammalian cell
Target Names	BZW2
Protein Names	Recommended name: Basic leucine zipper and W2 domain-containing protein 2
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.