



Recombinant Chicken Zinc finger Ran-binding domain-containing protein 2 (ZRANB2)

Product Code	CSB-EP736605CH
Abbreviation	ZRANB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLX5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MSTKNFRVSD GDWICPDKKC GNVNFARRTS CNRCGREKTT EAKMMKAGGT EIGKTLAEKS RGLFSANDWQ CKTCGNVNW RRSECNMCNT PKYAKLEERT GYGGGFNERE NVEYIEREES DGEYDEFGRK KKKYRGKPVG PASILKEVED KESEGEDEED EDGDLISKYKL DEDEDEDDAD LSKYNLDASE EEDTNKKKKS NRRSRSKSRS SHSRSSSRSS SHSSSRSRSR SHSRSSSSSR SRSRSSSREH SRSRGSKSRS SRSYRGSST PRKRSYSSSR SSSPERSKK RSRSRSSSSG DRKKRRSRSR SPERRRRSSS GSSHSGSRTS SKKK
Source	E.coli
Target Names	ZRANB2
Protein Names	Recommended name: Zinc finger Ran-binding domain-containing protein 2
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.