



# Recombinant Legionella pneumophila subsp. pneumophila Chromosomal replication initiator protein DnaA (dnaA)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP736664LDL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Abbreviation</b>      | dnaA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                     |
| <b>Uniprot No.</b>       | Q5ZZK8                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Immunogen Species</b> | Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Sequence</b>          | MSTTAWQKCL GLLQDEFSAQ QFNTWLRPLQ AYMDEQRLIL LAPNRFVVDW<br>VRKHFFSRIE ELIKQFSGDD IKASIEVGS KPVEAVDTPA ETIVTSSSTA<br>PLKSAPKKAV DYKSSHLNKK FVFDSEVEGN SNQLARAASM QVAERPGDAY<br>NPLFIYGGVG LGKTHLMHAI GNSILKNNPE AKVLYLHSER FVADMVKALQ<br>TNSINEFKRF YRSLNALLID DIQFFAGKDR SQEEFFHTFN ALLEGQQQII<br>LTSDRYPKEI EGMEERLKSR FGWGLTVAVE PPELETRVAI LISKAQSN<br>ELPYEVAFFI AKRIRSNVRE LEGALRRVIA NAHFTGKPIT IEFVHEALRD<br>LLALQDKLVT IENIQKTVAE YYKVKVADLL SKRRSRSIAR PRQMAMALSK<br>ELTNHSLPEI GDHFGGKDHT TVIHACRKVK ELIQDDSDFA EDYKNLMRIL SS |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | dnaA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Names</b>     | Recommended name: Chromosomal replication initiator protein DnaA                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Expression Region</b> | 1-452                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                            |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,                                                                                                                                                                                                                                                                                                                                                                                                                                                          |



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.