



Recombinant Mycoplasma hyopneumoniae

Elongation factor P (efp)

Product Code	CSB-BP736670MWK
Abbreviation	efp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q600M5
Product Type	Recombinant Protein
Immunogen Species	Mycoplasma hyopneumoniae (strain 232)
Purity	≥85% (SDS-PAGE)
Sequence	MINVNEFRPG ITFEFENEIY VVISAQHSKQ GRGQANVKAK VKNLRTGAIT IKTFSGGERV EKAHIEKISM SFLYNDGVSI VLMDDSTYEQ VAIENTKITW ELNFLTGEGIK VKLRKFNNEI LDIELPAKIE LKITSTFDAV RGNTTTNPTK RATLETGYEI DVPLFIKEGE SVIVSTEDGK YVSRA
Source	Baculovirus
Target Names	efp
Protein Names	Recommended name: Elongation factor P Short name= EF-P
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.