



Recombinant Methanocaldococcus jannaschii Translation initiation factor 6 (eif6)

Product Code	CSB-MP736703MRU
Abbreviation	eif6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60357
Product Type	Recombinant Protein
Immunogen Species	Methanocaldococcus jannaschii (strain ATCC 43067 / DSM 2661 / JAL-1 / JCM 10045 / NBRC 100440) (Methanococcus jannaschii)
Purity	≥85% (SDS-PAGE)
Sequence	MTMIIRKYFS GIPTIGVLAL TTEEITLLPI FLDKDDVNEV SEVLETKCLQ TNIGGSSLVG SLSVANKYGL LLPKIVEDEE LDRIKNFLKE NNLDLNV EII KSKNTALGNL ILTNDKGALI SPELKDFKKD IEDSLNVEVE IG TIAELPTV GSNAVVTNKG CLTHPLVEDD ELEFLKSLFK VEYIGKGTAN KGTTSSVGACI IANSKGAVVG GDTTGPELLI IEDALGLI
Source	Mammalian cell
Target Names	eif6
Protein Names	Recommended name: Translation initiation factor 6 Short name= aIF-6
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.