



# Recombinant *Methylococcus capsulatus* Glycogen synthase 2 (glgA2)

<b>Product Code</b>	CSB-MP736713MFM
<b>Abbreviation</b>	glgA2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q604D9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MHVTPELAPV AKVGGGLADV V FGLGRELEIR GNHVEIILPK YDCMRYDQIW GLQRTFDDLLW VPWYGGAIHC SVYFGFVHGR KCFEIEPHSQ DNFFNRGAVY GFHDDIFRFA FFSRAAMEFL WKAGKNPDII HCHDWQTALV PVYLYEIQP MGMRHPRVCF TIHNFKHQGV TGAQVLHASG LDRPEYYFHY DRLRDNHNPH AINLMKGGIV YANFVTTVSP RYAMEAKDQG QGFGLEPTLH IHMKEYGGV NGIDYDVWNP EIDPHIPVHF NVDTIEGKYA DKKALRDRL LADNEKPIVS FVGRLLDPQKG IELIRHALFY TLGQGGQFVL LGSSPDGAIN GYFWGLKRQF NDNPDCHLEI GYNEELAHLV YAGSDVMVVP SRFEPGLTQ LIAMRYGTIP VVREIGGLAD TVIDKDFSHR PLHERNGYVF RDYDERGLES ALGRAIACY QYPDHFRELM KNAMRYDYSW NHPGQDYLN I YHYIRDK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	glgA2
<b>Protein Names</b>	Recommended name: Glycogen synthase 2 EC= 2.4.1.21 Alternative name(s): Starch [bacterial glycogen] synthase 2
<b>Expression Region</b>	1-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.