



Recombinant Methylococcus capsulatus Type III pantothenate kinase (coaX)

Product Code	CSB-MP736748MFM
Abbreviation	coaX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q607I0
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MKLLIDIGNS RVKWARSEPG RYTFGGAFAT DASSLAEEACE SQMAALPRPE AVLVSNVAGM RAGRLLADWA SARWVVPVRF VRSEASALGV VNAYPNPAQL GVDRWLGLIG AWRTRRDACC IVSAGTALTI DLLDAEGRHR GGVIGPGLAM MRAALASLEA LGGAAESCDG FFCNRTGPAL SSGILHSAVG LIRESQARAA RLLATEPELL LTGGDAETLR PYLETPVRLV PHLVLEGLHA VAERMP
Source	Mammalian cell
Target Names	coaX
Protein Names	Recommended name: Type III pantothenate kinase EC= 2.7.1.33 Alternative name(s): PanK-III Pantothenic acid kinase
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.