



Recombinant Mouse Zinc finger and BTB domain-containing protein 17 (Zbtb17)

Product Code	CSB-MP736751MO
Abbreviation	Zbtb17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60821
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p>MDFPQHSQRV LEQLNQQRQL GLLCDCTFVV DGVDFKAHKA VLAACSEYFK MLFVDQKDVV HLDISNAAGL GQVLEFMYTA KLSLSPENV DVLAVASFLQ MQDIVTACHT LKSLAEPSST TGESADASAV EGGDKRAKDE KAAATMLSRL DQARGSSSTG PGRELKEERG GQAESASSGA EQTEKADAPR EPPPVELKPD PTSSMAAAEA EALSESSEQE MEVEPASKGE DGQEEEGAGP ATVKEEGMHL DN GEPPEENE ESAGTDSGQE LGMEGQNLRS GTYGDRTESK AYGSIHKCE DCGKEFTHTG NFKRHIRIHT GEKPFSCREC SKAFSDPAAC KAHEKTHSPL KPYGCEECKG SYRLISLLNL HKKRHSGEAR YRCGDCGKLF TTSGNLKRHQ LVHSGQKPYQ CDYCGRSFSD PTSKMRHLET HDTDKEHKCP HCDKKFNQVG NLKAHLKIHI ADGPLKCREC GKQFTTSGNL KRHLRIHSGE KPYVCTHCQR QFADPGALQR HVRIHTGEKP CQCVICGKAF TQASSLIAHV RQHTGEKPYV CERCGRFVQ SSQLANHIRH HDNIRPHKCS VCSKAFVNVG DLSKHIIIHT GEKPYLCDKC GRGFNRVDNL RSHVKT VHQG KAGIKILEPE EGGEVSVTV DDMVTLATEA LAATAVTQLT VVPVGA AVTA DETEVLKAEI SKAVKQVQEE DPNT HILYAC DSCGDKFLDA NSLAQH VRIH TAQALVMFQT DADFYQQYGP GSTWPAGQML QAGELVFRPR DGTEGQPTLA ESPPTAPDCL PPAE</p>
Source	Mammalian cell
Target Names	Zbtb17
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 17 Alternative name(s): LP-1 Polyomavirus late initiator promoter-binding protein Zinc finger protein 100 Short name= Zfp-100 Zinc finger protein 151 Zinc finge
Expression Region	1-794
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein

**Target Details**

This gene encodes a zinc finger protein involved in the regulation of c-myc. The symbol MIZ1 has also been associated with PIAS2 which is a different gene located on chromosome 18.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.