



Recombinant Mouse Inactive 2'-5'-oligoadenylate synthase 1B (Oas1b), partial

Product Code	CSB-MP736752MO
Abbreviation	Oas1b
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q60856
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MEQDLRSIPASKLDFIENHLPDTSFCADLREVIDALCALLKDRSFRGPVRRMR ASKGVKKGKGTTLKGRSDADLVVFLNNLTSFEDQLNQQGVLIKEIKKQLCEVQH ERRCGVKFEVHSLRSPNSRALSFKLSAPDLLKEVKFDVLPAYDLLDHLNLIKPP NQQFYANLISGRTPPGKEGKLSICFMGLRKYFLNCRPTKLKRLRLVTHWYQLC KEKLGDPPLPPQYALELLTVYAWEYGSRVTKFNATAQGFRTVLELVTKYKQLQIY WTVYYDFRHQEVSEYLHQQLKKDRPVILDPADPTRNIAGLNPKDWRRRLAGEA AAWLQYPCFKYRDGSSVCSWEVPTVEGVPMK
Source	Mammalian cell
Target Names	Oas1b
Protein Names	Recommended name: 2'-5'-oligoadenylate synthase 1B Short name= (2-5')oligo(A) synthase 1B Short name= 2-5A synthase 1B EC= 2.7.7.- Alternative name(s): 2'-5'-oligoadenylate synthase-like protein 1
Expression Region	1-351aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.