



Recombinant Methylococcus capsulatus Integration host factor subunit alpha (ihfA)

Product Code	CSB-MP736774MFM
Abbreviation	ihfA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60AY8
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MALTKADMVE RLFSELGLNK RDAKELVDQF FEVIKASLES GRSVKLSGFG NFDLRDKRQR PGRNPKTGEE IPITARRVVT FKAGQKLKAT VEALQEPELG KLDAERTGEV
Source	Mammalian cell
Target Names	ihfA
Protein Names	Recommended name: Integration host factor subunit alpha Short name= IHF-alpha
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.