



# Recombinant *Methylococcus capsulatus* ATP synthase subunit beta (atpD)

<b>Product Code</b>	CSB-EP736784MFM
<b>Abbreviation</b>	atpD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60CR4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSGKIVQII GAVVDAEFPR ESLPKIYDAL RVDEVGLVLE VQQQLGDGVV RCIAMGSTDG VRRGVTVTNT GAPISVPVGT ETLGRIMNVL GDPIDEKGPI GEKERWGIHR KAPSYQDQAA ANELLETGIK VIDLICPFAK GGVKVLFGGA GVGKTVNMME LIRNIAIEHS GYSVFAGVGE RTREGNDFYH EMKESNVLDK VSLVYQQMNE PPGNRLRVAL TGLTIAEKFR DEGRDVLLFI DNIYRYTLAG TEVSALLGRM PSAVGYQPNL AEEMGILQER ITSTKTGSIT SIQAVYVPAD DLTDPSPATT FAHLDATVVL SRQIAELGIY PAVDPLDSTS RQLDPLVIGH EHYETTRRVQ STLQRYKELR DIIAILGMDE LSENDKLIVS RARKIQRFLS QFFFVAEVFT GSPGKYVPLK DTIAGFKGIV DGDYDDVPEQ AFYMVGSIDE ALEKAKRLAA
<b>Source</b>	E.coli
<b>Target Names</b>	atpD
<b>Protein Names</b>	Recommended name: ATP synthase subunit beta EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit beta F-ATPase subunit beta
<b>Expression Region</b>	1-460
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.