



Recombinant Methylococcus capsulatus ATP synthase subunit beta (atpD)

Product Code	CSB-MP736784MFM
Abbreviation	atpD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60CR4
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MSSGKIVQII GAVVDAEFPR ESLPKIYDAL RVDEVGLVLE VQQQLGDGVV RCIAMGSTDG VRRGVTVTNT GAPISVPVGT ETLGRIMNVL GDPIDEKGPI GEKERWGIHR KAPSYQDQAA ANELLETGIK VIDLICPFAK GSKVGLFGGA GVGKTVNMME LIRNIAIEHS GYSVFAGVGE RTREGNDFYH EMKESNVLDK VSLVYQGMNE PPGNRLRVAL TGLTIAEKFR DEGRDVLLFI DNIYRYTLAG TEVSALLGRM PSAVGYPNL AEEMGILQER ITSTKTGSIT SIQAVYVPAD DLTDPSPATT FAHLDATVVL SRQIAELGIY PAVDPLDSTS RQLDPLVIGH EHYETTRRVQ STLQRYKELR DIIAILGMDE LSENDKLIVS RARKIQRFLS QPFVFAEVFT GSPGKYVPLK DTIAGFKGIV DGDYDDVPEQ AFYMVGSIDE ALEKAKRLAA
Source	Mammalian cell
Target Names	atpD
Protein Names	Recommended name: ATP synthase subunit beta EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit beta F-ATPase subunit beta
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.