



Recombinant Putative H/ACA ribonucleoprotein complex subunit 4 (CBG18316)

Product Code	CSB-EP736798CXX
Abbreviation	CBG18316
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60YA8
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MGKKDKRSLK EGDELAEAQQ KGSFQLPSSN ETAKLDASQW PLLLKNYDKL NVRTNHYTPH VEGVSPLKRD IKNYISSGFF NLDKPSNPSS HEVVSUWKRI LRCEKTGHSG TLDPKVSGCL IVCIDRTTRL AKSQQGAGKE YICIFKLHEE VEDERKVKQA LEKLTGALFQ RPPLISAVKR QLRIRTVYEN KFIEYDPAQQ MGIFNCICES GTYVRTICVH LGLILGCGGQ MQELRRNRSG ICDENENMVT MHDVLDAQYM LDTQKDESYM RHIVRPLEAL LTQHKRVVVK DSCVNAICYG AKILIPGILR YDDIEVGKE IVIMTTKGEA ICIAIAQMST STIASVDHGI VAKSKRVIME RDVYGRKWGL GPVASKKKQM VKDGLLDKFG KPNITTPKSW AKEYVQTDKV KKEQEDKEDE EEEEAPKKKS KKAACKEVSS SSDSE
Source	E.coli
Target Names	CBG18316
Protein Names	Recommended name: Putative H/ACA ribonucleoprotein complex subunit 4 EC= 5.4.99.-
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.