



Recombinant Skn-1 dependent zygotic transcript 1 protein (sdz-1)

Product Code	CSB-MP737059CXX
Abbreviation	sdz-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61QT3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MQFSIVYFLV FLIPSTQSTG FLSLKLTS DH DCLVHLEIQS MEYSETVRLL AYETKSLEIF TPEITNQLSI HFQILHHFSG VSLSDVSSQL FKLDNNGIWE PRVIDSDQVI LSIQSNFHCQ KGFYGPICER RSRSTSTIKT TTSIPTATSS PFHLPISVQ INNLIYAVL SSVVLLIIA NCFLCFCRPK PSKYIDYTFQ NVFPMDEFFP MEKTIGSPDT EYFDSSSRY Y STITLQSRV
Source	Mammalian cell
Target Names	sdz-1
Protein Names	Recommended name: Skn-1 dependent zygotic transcript 1 protein
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.