



# Recombinant Mouse Semaphorin-3B (Sema3b)

<b>Product Code</b>	CSB-MP737070MO
<b>Abbreviation</b>	Sema3b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q62177
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PNLPR LRLSFQELQA RHGVRTFRLE RTCCYEALLV DEERGRLFVG AENHVASLSL DNISKRAKKL AWPAPVEWRE ECNWAGKDIG TECMNFVKLL HTYNHHTLLA CGTGAFHPTC ALWRWATAGG THASTGPEKL EDGKGKTPYD PRHRAASVLV GEELYSGVTA DLMGRDFTIF RSLGQNP SLR TEPHDSRWLN EPKFVKVFWI PESENPD DDK IYFFFRESAV EAAPAMGRMS VSRVGQICRN DLGGQRSLVN KWTTFLKARL VCSVPGVEGD THFDQLQDVF LLSSRDRQTP LLYAVFSTSS GVFQGSACV YSMNDVRRAF LGPFAHKEGP THQWVSYQGR VPYPRPGMCP SKTFGTFSST KDFPDDVIQF ARNHPLMYP VLPMGGRPLF LQVGAGYTFT QIAADRVA A A DGHYDVLFIG TDVGTVLKVI SVPKGSRPNS EGLLLEELQV FEDSAAITSM QISSKRQQLY IASRSAVAQI ALHRCTALGR ACAECCLARD PYCAWDGSAC TRFQPTAKRR FRRQDIRNGD PSTLCSGDSS HSVLLEKKVL GVESGS AFLE CEPRSLQAHV QWTFQGAGEA AHTQVLA EER VERTARGLL RGLRRQDSGV YLCVAVEQGF SQPLRRLVLH VLSAAQAERL ARAEAAAPA PPGPKLWYRD FLQLVEPGGG GGANSLRMCR PQPGHHSVAA DSRRKGRNRR MHVSELRAER GPRSAAHW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Sema3b
<b>Protein Names</b>	Recommended name: Semaphorin-3B Alternative name(s): Semaphorin-A Short name= Sema A
<b>Expression Region</b>	26-748
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	The semaphorin/collapsin family of molecules plays a critical role in the guidance of growth cones during neuronal development. The secreted protein encoded by this gene family member is important in axonal guidance and has been shown to act as a tumor suppressor by inducing apoptosis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.