



Recombinant Ephrin-4 (efn-4)

Product Code	CSB-MP737084CXX
Abbreviation	efn-4
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q624D2
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	DEHTVHWN STNSMFRNRH PSIEVRLGDV VRFVCPDNEG RKNGEYLTVY EVSEFAMGEC ALESNSREVI KCGVDTNTEK IIRTHQLPIG ESREPPKNVA QFIRSVNPIP NGKEYQPGQT YYYITTSSGK PGGIGQQMYG LCVSKNMRLS MKVLSSQPTP SPSSKPARSR TDARRQEDFI TKSSAELMGG QEDESENDN AHLPRDLEI ATNPKFRRPS QFDQAAASAG VLDGQFLKVV QMAKEGKTGT FENDREVQKS AEKDAWDPIN RHYVADLMNS AYKNANDRVV YQREPDFLIH EEDISTNSLG YSSSS
Source	Mammalian cell
Target Names	efn-4
Protein Names	Recommended name: Ephrin-4
Expression Region	23-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>