



Recombinant Phodopus sungorus Thyrotropin subunit beta (TSHB)

Product Code	CSB-MP737085EUT
Abbreviation	TSHB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62590
Product Type	Recombinant Protein
Immunogen Species	Phodopus sungorus (Striped hairy-footed hamster) (Djungarian hamster)
Purity	>85% (SDS-PAGE)
Sequence	YCLTINTTIC AGYCMTRDIN GKLFLPKSAL SQDVCTYRDF TYRTVEIPGC PHHVAPYFSY PVAMSCKCGK CNTDYSDCIH EAVKTNCTK PQTFYLGGS V
Source	Mammalian cell
Target Names	TSHB
Protein Names	Recommended name: Thyrotropin subunit beta Alternative name(s): Thyroid-stimulating hormone subunit beta Short name= TSH-B Short name= TSH-beta Thyrotropin beta chain
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.